Paratec / Paratec Premium Instructions for use

Ready to use paraffin for laboratory use for the preparation of tissue or cellular specimens for microtomy. The specimens can then be evaluated by diagnostic techniques such as staining followed by study of images. Slee paraffin is DMSO-free and highly purified paraffin. The granulated material ensures an easy and convenient handling.

Only qualified and trained laboratory personnel should work with Slee paraffin.

For application, fill your beakers / containers with paraffin granules up to the desired level. After initial heating and melting down, further granules must be added to obtain the desired fill volume.

Only work with specimens that have been carefully and properly prepared. Only then the tissue can be processed / embedded / cut / stained and finally be diagnosed in the best possible way.

- In the tissue processor: All samples must be fixed, dehydrated and cleared prior to further processing.
- In the embedding station: All samples must be fixed, dehydrated and infiltrated with Paratec prior to further processing.

Store at room temperature (15 - 25 °C) in a dry area. Dispose of the used paraffin according to your laboratory instructions.

Expiration date: see packaging / label

Tips and tricks processing:

- Use an appropriate schedule
- Provide additional fixation
- Maintain reagent quality
- Use high-quality wax
- Avoid hazardous reagents

Tips and tricks embedding:

- Orientate specimen carefully
- Choose an appropriate mold and paraffin
- Handle specimen gently
- Avoid excessive heat
- Check temperature regularly
- Do not over-fill molds

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slee
solutions
for pathology

PARATEC

TISSUE PROCESSING PARAFFIN
DSMO-FREE

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
**PARATEC
PREMIUM**

TISSUE EMBEDDING PARAFFIN
FOR PRECISE EMBEDDING AND CUTTING
WITH POLYMER ADDITIVE
DSMO-FREE

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DESIGN &
MANUFACTURING
MADE IN GERMANY

Process of tissue preparation / Ablauf der Gewebepreparation

Features Paratec

Paratec was especially designed for **tissue processing** (fixation, dehydration, clearing and infiltration). Due to its molecular structure, it is particularly suitable for thorough and complete penetration of tissue samples.

Melting range 54 – 56 °C

Paratec is ideally suited for use in the following Slee devices:

- MTP carousel tissue processor
- MTM I and MTM II vacuum tissue processors

Features Paratec Premium

With polymer additive for **precise embedding and cutting**.

Paratec Premium was especially developed and is intended to be used in **tissue embedding**. Due to the long-chained polymer additive, **Paratec Premium** supports the stability of the specimen block especially well. Cutting will be easier and the sections will be significantly smoother when using **Paratec Premium** for embedding.

Melting range 54 – 57 °C

Paratec Premium is ideally suited for use in the following Slee devices:

- MPS P dispensing module
- MPS P I paraffin dispensing unit
- MPS P II modular paraffin embedding center

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Did you know?

Our Slee service team is there for you 24/7!

If your samples don't turn out the way you expected don't hesitate to contact us. We will be glad to help you. Fill out the contact form on our homepage:

www.slee.de/en/service/product-service/

We will get back to you asap and work it out together!

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SLEE medical GmbH · Am Neuberg 14