# Instructions for use

## Intended use

The Slee high quality blades are specifically designed for cutting thin sections of human tissue specimens, e. g. for cancer diagnosis.

## Hazard warnings

- Risk of injury. Handle with care!
- Blades are very sharp.  
They may only be used by trained laboratory personnel.
- Never touch the blade edge!
- Do not place blades at unsecured areas.
- Never position blades with the sharp edge pointing upwards.
- Store blades in a covered container.  
Use a container that has guides to hold the blades rigid.
- Never try to catch a falling blade.

## 1. Handling instructions

- Inspect the blades upon receipt. If the packaging is damaged, do not use the blades. Contact the manufacturer.
- If the blade is not being ejected from the dispenser, do not try to force the blade out. Contact the manufacturer.
- To push out the blade, pull the slider to the right (locking position).  
Then push it slowly all the way to the left, following the direction of the arrow.  
The blade will be pushed out of the left side of the dispenser.  
Take out the blade.

## 2. Disposal

- Insert the used blades (without residues or section waste) into the integrated waste box at the bottom of the dispenser. The opening is labeled with "USED BLADES".  
All blades will be stored safely and without danger to the laboratory staff.
- Make sure to insert the blades flat so all used blades will fit into the container.
- Please note that contaminated blades (from infectious or radioactive material) must be disposed of in a special container according to your lab requirements.

## 3. Technical data

### 3.1 Dimensions

- Item No.: 28407005 Low profile plasma LPS for soft specimens, plasma coated: 80 x 8 mm, thickness 0.254 mm, 35° angle
- Item No.: 28407004 Low profile plasma LPH for hard specimens, plasma coated: 80 x 8 mm, thickness 0.254 mm, 35° angle
- Item No.: 28407007 High profile plasma HPS for soft specimens, plasma coated: 80 x 14 mm, thickness 0.254 mm, 35° angle

### 3.2 Material

- Stainless steel

### 3.3 Transport and operating temperature / Storage humidity range

- Transport temperature: 10 °C to 40 °C
- Storage humidity: < 80 % non-condensing
- Operating temperature range: -40 °C to +40 °C

## 4. Warranty

SLEE medical GmbH guarantees that the product delivered has been subjected to a comprehensive quality control procedure, that the product is faultless and complies with all technical specifications and / or agreed characteristics warranted. Mistakes and errors which occur because of improper use cannot be accepted.

Questions? We are there for you!

## Tips and tricks cutting:

- Use high-quality blades
- Optimize knife tilt angle (clearance angle)
- Carefully trim the block
- Avoid freezing damage
- Use cold blocks
- Cut sections slowly
- Mind the specimen orientation

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PLASMA COATED  
PLASMABESCHICHTET

DESIGN &  
MANUFACTURING  
MADE IN GERMANY

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SLEE medical GmbH · Am Neuberg 14