

Tissue Embedding Paraffin

paratec premium

paratec premium consists of highly purified paraffin with the addition of polymer for precise embedding and sectioning. Melting range 54-57°C



<i>paratec premium</i>		
Description	Order #	Packing
High quality paraffin for precise tissue embedding	#30010004	8 x 1 kg
	#30010006	4 x 2.5 kg
	#30010007	10 kg

SLEE medical GmbH

Lise-Meitner-Strasse 11
55129 Mainz · Germany

Tel. +49 6131 95871-0
Fax +49 6131 95871-722


mail@slee.de
www.slee.de

Brochure_consumables_en_2018-10-003.00

SLEE

 medical

Consumables



- High quality plasma blades
- Slides
- Cover slips
- Embedding Paraffin


DESIGN AND MANUFACTURING
MADE IN GERMANY




High Quality plasma blades

LPS, LPH, HPS


The new SLEE blades for microtomes feature a special plasma technology, improving the sharpness of the blades as well as the durability of the blades considerably.



LPS		
Description	Order #	Packing
Low profile plasma blade LPS for soft specimens	#28407005	50 pc/box



LPH		
Description	Order #	Packing
Low profile plasma blade LPH for hard specimens	#28407004	50 pc/box



HPS		
Description	Order #	Packing
High profile plasma blade HPS for soft specimens	#28407007	50 pc/box

High Quality slides / cover slips

The new SLEE StarFrost® cover slips in 3 different sizes, made for automated cover-slipping and the new SLEE StarFrost® slides in different styles and colors perfectly matching the SLEE printtec s slide printer.



Slides		
Description	Order #	Packing
StarFrost® Slides 26 x 76 x 1.0 mm, ground edges 90°, color frosted marking area about 20 mm (available in white, orange, blue, green, yellow, pink)	On request	50 pc/box
StarFrost® Slides 26 x 76 x 1.0 mm adhesive, ground edges 90°, color frosted marking area about 20 mm (available in white, orange, blue, green, yellow, pink)	On request	50 pc/box
StarFrost® Slides 26 x 76 x 1.0 mm advanced adhesive, ground edges 90°, color frosted marking area about 20 mm (available in white, orange, blue, green, yellow, pink)	On request	50 pc/box



Cover slips		
Description	Order #	Packing
StarFrost® Cover Slips 24 x 50 mm	#30001051	100 pc/box
StarFrost® Cover Slips 24 x 55 mm	#30001052	100 pc/box
StarFrost® Cover Slips 24 x 60 mm	#30001053	100 pc/box