

Product specification

Ergotec · Worktable

Item No.: 10090056
Standard version with height adjustment
1,800 x 800 mm



Specifications

Range of height adjustment	660 up to 1,100 mm
Programmable height levels	4 memory settings
Maximum load	400 kg
Thickness of table top	30 mm
Color of table top	selectable
Maximum dimensions (possible customization)	2,000 x 1,000 mm
Cooling plate dimensions	300 x 320 mm
Cooling plate temperature	+15 °C down to -15 °C
Water bath dimensions	260 x 214 x 60 mm
Water bath temperature	RT up to 85 °C

Standard components

Ergotec - standard version with height adjustment 1,800 x 800 mm

Manual

Options*

	Item No.
Integrated water bath	10120002
Integrated cooling plate	10120003
Cable channel*	10120004
Workplace with LED illumination*	10120005
Undertable drawer*	10120006
Ergonomic cutout	on demand

*Contact us for further information.

Configuration variants*

Dimension of table (possible width: 80, 100, 120, 140, 160, 180, 200 cm / possible depth: 80, 90, 100 cm)

Color of table top

Edge shape of table top (straight edge / bevel edge)

Corner radius of table top (1 cm, 2 cm, 3 cm, 4 cm)

Adjustable height - motorized / static height

Set-up of sockets and network in cable channel

Color and shape of cable channel

*Contact us for further information. Due to numerous configuration variants, we offer a product configuration sheet, available for download at www.slee.de.

Optional accessories

	Item No.
Ergosit laboratory chair	10157010

*Contact us for further information.

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
E: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)