


# Product specification

## Disposable blades · plasma coated

Item No.: 28407005      Low profile plasma LPS  
 Item No.: 28407004      Low profile plasma LPH  
 Item No.: 28407007      High profile plasma HPS



### Description

Disposable blade - Item No.	28407005	28407004	28407007
Profile	low (LPS)	low (LPH)	high (HPS)
Dimensions (W x D x H)	80 x 8 x 0.254 mm	80 x 8 x 0.254 mm	80 x 14 x 0.254 mm
Blade angle	35°	35°	35°
Blade hardness	600 Hv	600 Hv	600 Hv
Edge width	20 µm	20 µm	20 µm
Allowable edge blunting	± 1°	± 1°	± 1°
Material	stainless steel	stainless steel	stainless steel
Blades per dispenser	50 blades	50 blades	50 blades
Application	for soft tissues	for hard tissues	for hard tissues

### Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14  
 55268 Nieder-Olm · Germany  
[www.slee.de](http://www.slee.de)

T: +49 (0) 6136 76997-0  
 E: [mail@slee.de](mailto:mail@slee.de)



[www.slee.de](http://www.slee.de)



[YouTube](#)



[LinkedIn](#)



[Instagram](#)