




CONSUMABLES

BLADES WITH PLASMA COATING
BLADES WITH PTFE COATING
PARAFFIN
SLIDES / COVERSLIPS
CRYOGLUE EMBEDDING MEDIUM
ACCESSORY FOR CRYOSTATS

DESIGN &
MANUFACTURING
MADE IN GERMANY

Consumables

All Slee consumables are subject to the highest quality requirements. Our product range includes blades with plasma coating, blades with PTFE coating, tissue embedding paraffin, slides and coverslips as well as the CryoGlue embedding medium and other consumables.



High quality microtome blades

LPS, LPH, HPS, PTFE coated


The Slee blades for microtomes feature a special plasma technology, that considerably improves the sharpness as well as the durability of the blades. In addition, further blades with innovative PTFE coating are available. All Slee microtome blades are made of surgical steel.

Disposable blades (plasma coated)

Item No.

Low profile plasma LPS
for soft specimens, 50 pcs. / box

28407005



Disposable blades (PTFE coated)

Item No.

Low profile suitable for all types
of specimens, 50 pcs. / box


28407000



Low profile plasma LPH

for hard specimens, 50 pcs. / box

28407004



High profile

suitable for all types
of specimens, 50 pcs. / box


28408000



High profile plasma HPS

for soft specimens, 50 pcs. / box

28407007



Paratec Premium, Paratec Tissue embedding paraffin

Paratec Premium consists of highly purified paraffin. With polymer additive for precise embedding and cutting, melting range 54 - 57 °C. SLEE paraffin is DMSO-free.

Paratec Premium	Item No.	Paratec	Item No.
Premium paraffin		Paraffin	
1 kg	30010004	without polymer additive	
2.5 kg	30010006	melting range 54 - 56 °C	
10 kg	30010008	20 kg	30010002



Slides / coverslips and accessory in premium quality

Perfect accessory for the Slee **Printtec S** slide printer / **MCS I & MCS II** coverslippers, **MEV / MEV +** and **MNT** cryostats.

Slides

Item No.

StarFrost® slides

26 x 76 x 1.0 mm,
ground edges 90°,
clipped slide corners 45° angle,
color frosted marking area about 20 mm
(available in white, orange, blue, green,
yellow, pink)
50 pcs. / box

on demand

StarFrost® slides

26 x 76 x 1.0 mm, adhesive,
ground edges 90°,
clipped slide corners 45° angle,
color frosted marking area about 20 mm
(available in white, orange, blue, green,
yellow, pink)
50 pcs. / box

on demand

StarFrost® slides

26 x 76 x 1.0 mm, advanced adhesive,
ground edges 90°,
clipped slide corners 45° angle,
color frosted marking area about 20 mm
(available in white, orange, blue, green,
yellow, pink)
50 pcs. / box

on demand

Coverslips

Item No.

StarFrost® coverslips

24 x 50 mm, 100 pcs. / box

30001051

StarFrost® coverslips

24 x 55 mm, 100 pcs. / box

30001052

StarFrost® coverslips

24 x 60 mm, 100 pcs. / box

30001053

Accessory for cryostats

Item No.

CryoGlue embedding medium (4 x 125 ml)

30001100

Cryostat low temperature oil, 50 ml

30001011

Section waste cloth (set with 10 pcs.)

30000015

Cork plates (Ø 20 mm; 100 pcs.)

30001001

Highlights

- Slee StarFrost® slides in different styles and colors
- Slee StarFrost® coverslips in 3 different sizes, made for automatic coverslipping

Service

We are committed to provide world-class support to our global customers at all times. We strive to advise you professionally before purchase and to support you during installation and initial setup. In addition, it goes without saying that we will always be available as a reliable partner during the longtime use of our products.

You can contact us directly with any kind of questions. In order to process your inquiries in time, please always include the following information:

- User name
- Address
- Phone number
- E-mail
- Type of device
- Serial number of device
- Date of purchase
- Detailed case description

Service Hotline: +49 (0) 6136 76997-0

Service E-mail: service@slee.de

Trainings

Service and application training

Well-trained staff is crucial for the implementation of regular service and maintenance of high quality.

We offer regular service and application training for our worldwide distribution partners.

Visit www.slee.de to access the training schedule and to download the registration form.

Training requests

Tel.: +49 (0) 6136 76997-0

E-mail: mail@slee.de

Notes



SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0

E: mail@slee.de



www.slee.de



YouTube



LinkedIn



Instagram

affin
ectic