


slee solutions
for
pathology

WATER BATHS
& SLIDE WARMERS
Slidetec Water / Heat
Slidetec Heat
Slidetec Water
Slidetec Water S

DESIGN &
MANUFACTURING
MADE IN GERMANY

Overview

Water baths & slide warmers



Slidetec Water / Heat

Water bath & slide warmer

Specifications

Water bath volume	2 l
Capacity slide warmer	30 slides
Temperature range	RT up to 60 °C (water bath) RT up to 80 °C (slide warmer)
Electrical parameters	115 V, 60 Hz 230 V, 50/60 Hz
Consumption	600 W
Dimensions (W x D x H)	370 x 560 x 230 mm
Weight	18 kg

Item No.

Slidetec Water / Heat 230 V, 50/60 Hz 13000500

Slidetec Water / Heat 115 V, 60 Hz 13000501

Slidetec Heat

Slide warmer

Specifications

Capacity slide warmer	30 slides
Temperature range	RT up to 80 °C
Electrical parameters	115 V, 60 Hz 230 V, 50/60 Hz
Consumption	300 W
Dimensions (W x D x H)	350 x 370 x 230 mm
Weight	15 kg

Item No.

Slidetec Heat 230 V, 50/60 Hz 13000510

Slidetec Heat 115 V, 60 Hz 13000511



Slidetec Water

Water bath

Specifications

Water bath volume	2 l
Temperature range	RT up to 60 °C
Electrical parameters	115 V, 60 Hz 230 V, 50/60 Hz
Consumption	300 W
Dimensions (W x D x H)	350 x 370 x 230 mm
Weight	10 kg

Item No.

Slidetec Water 230 V, 50/60 Hz 13000520

Slidetec Water 115 V, 60 Hz 13000521

Optional accessories

Item No.

Lid for Slidetec Water or Water / Heat to cover the glass tray 13000550

Slidetec Water S

Water bath

Specifications

Water bath volume	1.6 l
Temperature range	RT up to 60 °C
Electrical parameters	115 V, 60 Hz 230 V, 50/60 Hz
Consumption	300 W
Dimensions (W x D x H)	280 x 280 x 100 mm
Weight	2 kg

Item No.

Slidetec Water S 230 V, 50/60 Hz 13000530

Slidetec Water S 115 V, 60 Hz 13000531

Optional accessories

Item No.

Lid for Slidetec Water S 13000551



SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0

E: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)