

# Product specification

## MPS C · Cooling module

Item No.: 10175000 230 V, 50/60 Hz  
 Item No.: 10177000 115 V, 60 Hz



## Specifications

Nominal power:	230 V AC + 10%
Nominal frequency	50 / 60 Hz
Power:	150 W
Protective class:	(1) I
Fuses	2 x T 6.3 A
Pollution degree:	(1) 2
Overvoltage protection category	II
Max. Heat emission	150 J/s
Operating temperature range	+10 to +35 °C
Operating humidity	max. rel. 80 % non-condensing
Storage temperature range	+5 to +55 °C
Storage humidity	max. rel. 80 % non-condensing
Dimensions [W x D x H]:	400 mm x 575 mm x 295 mm
Weight [w/o accessories]	24 kg
Temperature range	down to -15 °C
Cooling gas:	R 134 a
Work area	anodized aluminum 270 x 370 mm 80 cassettes

Hints for installation < 2000 m NN, 15 cm distance from backside to wall.

## Optional accessories

### Item No.

Transparent hood for cooling area MPS C

10175100

## Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14  
55268 Nieder-Olm · Germany  
[www.slee.de](http://www.slee.de)

T: +49 (0) 6136 76997-0  
M: [mail@slee.de](mailto:mail@slee.de)



[www.slee.de](http://www.slee.de)



[YouTube](#)



[LinkedIn](#)



[Instagram](#)