

## Product specification

### Slidetec Heat · Slide warmer

Item No.: 13000510      230 V, 50/60 Hz  
Item No.: 13000511      115 V, 60 Hz



### Specifications

Temperature range	up to 80 °C
Slide warmer capacity	30 microslides
Control	microprocessor, programming of temperature, working day, working time, time and date
Dimensions [W x D x H]	350 x 370 x 230 mm
Weight (w/o accessories)	15 kg
Electrical parameters	Item No.: 13000510    230 V; 50/60 Hz, Pmax: 300 W Item No.: 13000511    115 V; 50/60 Hz, Pmax: 300 W

### Standard components

Slide warmer (stage)
Operation manual
Mains cable

## Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14  
55268 Nieder-Olm · Germany  
[www.slee.de](http://www.slee.de)

T: +49 (0) 6136 76997-0  
M: [mail@slee.de](mailto:mail@slee.de)



[www.slee.de](http://www.slee.de)



[YouTube](#)



[LinkedIn](#)



[Instagram](#)