


Product specification

Slidetec Water / Heat · Water bath & slide warmer

Item No.: 13000500 230 V, 50/60 Hz
Item No.: 13000501 115 V, 60 Hz



Specifications

Temperature range	RT up to 60 °C, up to 80 °C with lid for water bath RT up to 80 °C for slide warmer
Water bath capacity	2 liters
Water bath illumination	LED (white / blue)
Slide warmer capacity	30 slides
Control	microprocessor, programming of temperature, working day, working time, time and date
Dimensions [W x D x H]	560 x 370 x 230 mm
Weight (w/o accessories)	18 kg
Electrical parameters	Item No.: 13000500 230 V; 50/60 Hz, Pmax: 600 W Item No.: 13000501 115 V; 50/60 Hz, Pmax: 600 W

Standard components

Histology water bath (glass, LED illumination)

Slide warmer (stage)

Operation manual

Mains cable

Optional accessories

Item No.

Lid for Slidetec Water or Water / Heat to cover the glass tray

13000550

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
M: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)