

Product specification

CS II · Cytocentrifuge

Item No.: 15500100	230 V, 50/60 Hz Open rotor
Item No.: 15500101	230 V, 50/60 Hz Closed rotor
Item No.: 15501100	230 V, 50/60 Hz Open rotor, with cooling
Item No.: 15501101	230 V, 50/60 Hz Closed rotor, with cooling
Item No.: on demand	115 V, 60 Hz



Specifications

Power supply:	200 – 240 V 50 – 60 Hz
Power Consumption:	400 VA
Chamber volume:	Cellfunnel Single 0.5 ml per chamber Double 0.5 ml per chamber ECOfunnel 5 ml
Rotor Slide Capacity:	Cellcliprotor up to 12 microslides
Max. RPM (speed):	Cellcliprotor 2.000 min ⁻¹ Centrifuge with applicable rotor 6.000 min
Running time:	1 – 99 minutes, ∞ continuous run, short cycle mode (impulse key)
Safety:	Imbalance switch-off Automatic rotor recognition Lid locking and holding
Dimensions [W x D x H]:	395 x 520 x 346 mm
Weight [w/o accessories]:	31 kg

Standard components

	CS II Open rotor	CS II Closed rotor
Cellclip rotor, open (12 positions)	✓	
Cellclip rotor, closed (12 positions)		✓
Cellclip (stainless steel, 12 pcs.)	✓	✓
Cytofunnel, reusable (12 pcs.)	✓	
Filter cards for cytofunnel, reusable Item No. 39910001 standard (200 pcs.)	✓	✓
Cellslides Item No. 39910003 (100 pcs.)	✓	✓
Mains cable	✓	✓
Manual	✓	✓

Optional accessories

	Quantity	Item No.
Cellclip	1 pcs.	39910031
Cytofunnel single, reusable		
Filter cards (only for cytofunnel single 39910000) up to 0.5 ml	200 pcs.	39910001
Filter cards (only for cytofunnel single 39910000) thin, up to 0.4 ml	200 pcs.	39910006
Cellslides single, uncoated	100 pcs.	39910002
Cellslides single, coated	100 pcs.	39910003
Cytofunnel double, reusable		
Cytofunnel double	12 pcs.	39910010
Filter cards for cytofunnel double up to 2 x 0.5 ml	500 pcs.	39910011
Cellslides double, uncoated (only for cytofunnel double, reusable 39910010)	100 pcs.	39910012
Cytofunnel ECOfunnel®, reusable		
Cytofunnel ECOfunnel®	12 pcs.	39910020
ECOseal seal for ECOfunnel®	100 pcs.	39910021
ECOSlides, coated	100 pcs.	39910022
Rotor (closed) CS centrifuge, removable, 12 positions	1 pcs.	39910100
Rotor (open) CS centrifuge, removable, 12 positions	1 pcs.	39910050
EASY rotor (closed) 12 positions, without stainless steel clips, for EASY / double cytofunnel	1 pcs.	15500600
EASY rotor (open) 12 positions, without stainless steel clips, for EASY / double cytofunnel	1 pcs.	15500601
Swing-out rotor 8-point for 15 ml	1 pcs.	39910042

Optional accessories

	Quantity	Item No.
Swing-out rotor 6-point for 50 ml	1 pcs.	39910043
Swing-out rotor 6-point Adapter for 15 ml (only for 39910043)	1 pcs.	39910044
Swing-out rotor 6-point for 50 ml	1 pcs.	39910043
Swing-out rotor 6-point Adapter for 15 ml (only for 39910043)	1 pcs.	39910044
EASY cytofunnel single, disposable, with white filter card up to 0.5 ml	40 pcs.	39910090
EASY cytofunnel single, disposable, with brown filter card up to 0.5 ml		on demand
EASY cytofunnel double, disposable, with white filter card up to 2 x 0.5 ml	40 pcs.	39910091
EASY cytofunnel double, disposable, with brown filter card up to 2 x 0.5 ml		on demand
Cytofunnel single, disposable, up to 0.5 ml	100 pcs.	39910004
Cytofunnel single, disposable, up to 0.5 ml	500 pcs.	39910005
Cytofunnel single, disposable, uncoated	100 pcs.	39910002
Cytofunnel single, disposable, coated	100 pcs.	39910003
Cytofunnel double, disposable, sample chamber with 2 shifted holes, up to 2 x 0.5 ml	100 pcs.	39910014
Cytofunnel double, disposable, sample chamber with 2 shifted holes, up to 2 x 0.5 ml	500 pcs.	39910015
Cellslides double, uncoated (only for cytofunnel double 39910015)	100 pcs.	39910016
Cellslides double, coated (only for cytofunnel double 39910015)		on demand
ECOfunnel®, disposable, up to 5 ml	100 pcs.	39910023
ECOfunnel®, disposable, up to 5 ml	500 pcs.	39910024
ECOfunnel®, disposable, ECOslides, coated	100 pcs.	39910022
ECOfunnel®, disposable, ECOslides, uncoated		on demand

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
M: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)