

## Product specification

### Slidetec Water S · Water bath

Item No.: 13000530      230 V, 50/60 Hz  
Item No.: 13000531      115 V, 60 Hz



### Specifications

Exterior dimensions (diameter / height)	280 mm / 100 mm
Interior dimensions (diameter / height)	200 mm / 60 mm
Usable bath height	50 mm
max. filling quantity (in liters)	1.6 l
Temperature range	approx. 5 °C above ambient to +80 °C
Temperature control	capillary thermostat
Temperature constancy, temporal	approx. ± 0.5
Power	220 - 240 V / 50-60 Hz
	100 - 115 V / 50-60 Hz
Protection type / protection class	IP20/I
Ambient conditions	use only indoors (not in potentially explosive surroundings)
Height above sea level	up to 2.000 m above sea level
Ambient temperature	+ 10 °C to + 40 °C
Humidity	max. 80 % relative humidity up to 31 °C
Weight	2 kg

## Standard components

Mains cable

---

Holder for thermostat

---

Thermostat

---

## Optional accessories

**Item No.**

Lid for Slidetec Water S

---

13000551

---

## Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14  
55268 Nieder-Olm · Germany  
[www.slee.de](http://www.slee.de)

T: +49 (0) 6136 76997-0  
M: [mail@slee.de](mailto:mail@slee.de)



[www.slee.de](http://www.slee.de)



[YouTube](#)



[LinkedIn](#)



[Instagram](#)