

## Product specification

### Disposable blades · PTFE coated

Item No.: 28407000      Low profile  
Item No.: 28408000      High profile



### Description

Disposable blade - Item No.	28407000
Profile	Low
Dimensions	80 x 8 mm
Blade Angle	35°
Material	Stainless steel
Blades per dispenser	50 blades
Blades per case	50 blades
Application	for all tissue types

### Description

Disposable blade - Item No.	28408000
Profile	High
Dimensions	77 x 14 mm
Blade Angle	35°
Material	Stainless steel
Blades per dispenser	50 blades
Blades per case	50 blades
Application	for all tissue types

## Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14  
55268 Nieder-Olm · Germany  
[www.slee.de](http://www.slee.de)

T: +49 (0) 6136 76997-0  
M: [mail@slee.de](mailto:mail@slee.de)



[www.slee.de](http://www.slee.de)



[YouTube](#)



[LinkedIn](#)



[Instagram](#)